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Il fabbricante: **FARO S.p.A.**
Via Faro, 15
20060 Ornago (Milano) ITALY

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CONTEXT

- Symbols used in the manual page 21
- Safety norms page 22-23-24-25
- Features page 26
- Description of the parts page 27
- Installation page 28-29
- Instruction for use page 30
- Maintenance page 31-32
- Troubleshooting page 33
- Technical specifications page 34-35

SYMBOLS USED IN THE MANUAL

Danger

Paragraphs marked with this symbol contain instructions to be followed carefully to avoid damage to the device, the operator and patient.

Warning

These instructions advise that great care must be taken to avoid situations that could damage the system.

Read carefully the instructions before use.

Forbidden

This icon points out what to prevent doing to not damage the device.

Suggestions

Information for most effective device use.

Dispose of the system in accordance with regulations for the separate collection of electrical material.

SAFETY NORMS

- Follow all instructions shown in the following manual.
- The motor must be powered by a FARO PN 617300 dedicated card; The wrong installation can influence the correspondence to the safety norms; for installation see card manual.
- Installation of the system must only be performed in a dental studio by specialized staff.

- **In order to avoid cross contamination between patients, insert or disconnect the dental studio micro-motor from the handpiece, avoiding contact with contaminated hands and gloves.**
- **The device is not protected by possible liquid penetration.**
- The device is not supplied as sterilized, therefore it must be cleaned before use (See page 28).

- **The device cannot be sterilized.**
- **According to EN 60825-1, don't point the led light directly into eyes, since the device belongs to class II.**

- The dental studio micro-motor “**Modo Infinity**” does not require lubrication.
- The device is connected to a “dental handpiece” conform to **UNI EN 23964**.

- **During final installation check the dielectric rigidity and the leakage currents along with EN 60601-1.**
- Do not attempt to repair the electrical parts. Do not intervene directly either on the dental studio micro-motor or on the corresponding electronic circuit.
- Carry out only the operations explained in this manual. In any other case contact the customer service.
- The device in its original packaging can be transported or kept in storage for a period of 15 weeks if the following environmental conditions are observed:
 - ⇒ **Ambient temperature from -20°C to $+70^{\circ}\text{C}$**
 - ⇒ **Relative humidity from 10% to 90%**
 - ⇒ **Atmospheric pressure from 500 to 1060 mBar.**

- The device is not protected from the effects of defibrillator electric shocks.

REQUIREMENTS FOR ELECTROMAGNETIC COMPATIBILITY

This medical device requires particular precautions regarding electromagnetic compatibility and must be installed and used in accordance with the information supplied with the accompanying documents.

| Guide and manufacturer's declaration – Electromagnetic Emissions | | |
|---|------------------|--|
| The dental studio micro-motor “Modo Infinity” can be used in the electromagnetic environment specified as follows. The customer or user must ensure that the device is used in that environment | | |
| Emissions test | Compliance | Guide to the electromagnetic environment |
| Irradiated Emissions CISPR11 | Class B , Group1 | The dental studio micro-motor “modo infinity” uses RF energy only for its internal operation. The RF emissions are therefore very low and do not cause any interference when near any other electronic device. The dental studio micro-motor “modo infinity” is suitable for use in dental studios and environments connected directly to the public mains supplying power to those environments. |
| Conducted emissions CISPR11 | Class B , Group1 | |
| Harmonic IEC/EN 61000-3-2 | Class A | |
| Harmonic IEC/EN 61000-3-3 | Conform | |

| Guide and manufacturer's declaration – Electromagnetic Immunity | | |
|---|---|---|
| The dental studio micro-motor “Modo Infinity” can be used in the electromagnetic environment specified as follows. The customer or user must ensure that the device is used in that environment | | |
| Emissions test | Compliance | Guide to the electromagnetic environment |
| Electrostatic discharges (ESD) IEC/EN61000-4-2 | $\pm 6\text{kV}$ by contact $\pm 8\text{kV}$ in the air | Floors should be in wood, cement or ceramic tiles. If floors are covered in synthetic material, the relative humidity must be a maximum of 30%. |
| Fast transients/bursts IEC/EN61000-4-4 | $\pm 2\text{kV}$ | The power supply should be that of a commercial or hospital environment. |
| Surge IEC/EN61000-4-5 | $\pm 1\text{kV}$ differential modus $\pm 2\text{kV}$ traditional modus | The power supply should be that of a commercial or hospital environment. |
| Voltage dips, brief interruptions and voltage variations IEC/EN61000-4-11 | $< 5\%U_T$ per 0.5 cycle of cycle $40\%U_T$ per 0.5 cycle of cycle $70\%U_T$ per 25 cycle of cycle $< 5\%U_T$ per 25 cycle of cycle | The power supply should be that of a commercial or hospital environment. If the user of the dental studio micro-motor needs the device to be used continuously, an uninterruptible power supply is recommended. |
| Magnet field IEC/EN61000-4-8 | 3A/m | The magnetic field should be that of a commercial or hospital environment. |
| Conducted immunity IEC/EN61000-4-6 | 3V_{rms} 150kHz to 80MHz | -- |
| Irradiated immunity IEC/EN61000-4-3 | 3V/m 80MHz to 2.5GHz | -- |
| Note U_T is the power supply voltage value | | |

FEATURES


Dental studio micro-motor versions “Modo Infinity”:

Intramatic led 40IL – Intramatic led 40ILU


- The dental studio micro-motor series “**Modo Infinity**” supply kinetic energy to a handpiece (not supplied) to be connected to its end. The device allows the passage of water and air through the cooling spray of the cutter, and gives light to the handpieces provided with optic fiber.
- The **dedicated** cord allows to supply the “**Modo Infinity**” version, model 40IL (intramatic-led)
- The **unified** cord micro-motor/turbine allows to supply power to the “**Modo Infinity**”, model 40ILU (unified-intramatic-led)

DESCRIPTION OF THE PARTS

Dimensions

- 1 Coupling to handpiece
- 2 O-ring
- 3 Air
- 4 Water

INSTALLATION

Check that the package contains the following components:

- Micro-motor
- 3 O-rings
- Instruction manual

- **THE DEVICE CAN BE INSTALLED ON THE USER'S CORD.**

- **"40IL" model**

Insert the micro-motor into the **dedicated cord** junction block.

Screw the ring nut of the feeding cord to the micro-motor, by rotating it clockwise.


- **"40ILU" model**

Unscrew the turbine ring nut from the **FARO unified cord** and insert the micro-motor into the junction block of the cord, by screwing it to the ring nut.

INSTALLATION

- **"40ILU"** model

Unscrew the ring nut turbine "A" from the FARO unified cord, and insert the micro-motor into the junction block of the cord.

INSTRUCTIONS FOR USE

After having connected the micro-motor to the cord, ensure its good functioning as follows:


- In the **40IL-40ILU** versions it should come out of the holes in correspondence with the O-rings, **AIR** from hole “**A**” and **WATER** from hole “**B**” (fig. 1). Led should go on. After this, engage the handpiece and check that there are no water or air leaks on the cord coupling.

IMPORTANT

Do not insert the handpiece when the micro-motor rotates.

fig. 1

MAINTENANCE

- **During installation the device must be cleaned.**
- **The device must be cleaned after every surgical intervention and before every use.**
- **Before cleaning the device, disconnect from the feeding cord.**
- **The device is not protected by liquid penetration.**
- **Clean the micro-motor, including the cord, with a cloth soaked in a disinfectant or in a non abrasive detergent.**

THE DEVICE CANNOT BE STERILIZED.

TO MAINTAIN QUALITY AND CORRESPONDENCE TO STANDARDS USE ONLY FARO SPARE PARTS.

Lubrication

It is suggested that in the “**Modo 40II and 40ILU**” versions the three “**A**” gaskets should be regularly lubricated with vaseline grease to limit wearing and to guarantee a perfect watertight and airtight capacity (fig. 2). Should they be worn-out, they must be replaced.

IMPORTANT

Do not lubricate any other mechanical parts of the micro-motor.


fig. 2

GUIDE TO THE MOST COMMON PROBLEMS – TROUBLESHOOTING

- The micro-motor has water or air leaks.

- The micro-motor does not rotate.

- The led does not go on.

- Air and water leaks from the cord junction

- Check the O-rings located on the handpiece junction. Should they be worn-out, replace them (see page 32 fig.2).
- If the problem persists, consult FARO technical assistance.
- Verify that the cord is correctly engaged.
- If the problem persists, consult FARO technical assistance.
- Consult FARO technical assistance.
- Verify the motor/cord coupling.
- If the problem persists, consult FARO technical assistance.

TECHNICAL SPECIFICATIONS

"Intramatic 40IL / 40ILU" Version

| | | |
|--|---|------------------------------------|
| Rotation speed | : | 1.500 to 40.000 rev/min. |
| Rotation | : | clockwise and anticlockwise |
| Led for light | : | 40lm @ 350mA |
| Max. absorbed power (with dedicated circuit) | : | 2 A |
| Given power | : | 45 W |
| Shaft max. torque | : | 10,8 Nm |
| Air consumption at the handpiece at 2,5 bar | : | 4,2 l/min |
| Water consumption of the spray at bar 2 bar | : | 160 ml/min |
| Consumption of motor cooling air at 2,5 bar | : | 12 l/min |
| Weight | : | 75 g |

Service with intermittent S3 load (CEI 34-1)

2 A 60 sec. ON - 10 min OFF

1,5 A 90 sec. ON - 12 min OFF


1 A 120 sec. ON - 6 min OFF

Class II device

Applied part type BF

Labeling conform to EN 60825-1

GUARANTEE CERTIFICATE

FARO offers the final customer a **12 month** guarantee starting from the date of purchase.

Repairs under guarantee must be performed at FARO; expenses and transport risks are at the risk of the purchaser.

Repair under guarantee is considered valid only when:

- **all sections of the certificate have been filled in and sent in advance to FARO by Fax (039.6010540).**

The guarantee covers faults due to the bad quality of the material or manufacturing defects; in the case of valid claims, the guarantee covers free repair or replacement. **Claims for damages and/or interest are excluded.**

The guarantee is not considered valid, at the sole discretion of FARO, if the fault is due to tampering, damage, incorrect use, improper maintenance and normal wear and tear.

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